

# BI@ BS5+

## Handheld Digital Medical Scope Series

### Product Features:


1. Specially designed for the department of ENT. It can be used to examine eardrum, nasal cavity and throat with different camera probe.
2. Avoid high cost of the traditional ENT examination equipment, such as light source, CCD, Monitor etc.
3. Light weight and portable, suitable for outpatient department, emergency room and medical ward.
4. High resolution image. The SD card can store multiple freeze images and videos, and can be connect to computer to print patient reports.
5. 3.5 inch high definition LCD screen, watching and taking photos at the same time.
6. It can be connected to large size monitor, synchronized to the big screen

### Product Appearance:



### Product specifications:

- Resolution: 720 x 480
- Light source: Natural white (LED)
- LCD screen: 3.5 inch high-definition LCD
- Image format: JPEG, BMP
- Media: Micro SD flash card
- Power supply: rechargeable lithium battery
- External power: 110~220 VAC, 50/60 Hz



**Wuxi Biomedical Technology Co., Ltd.**

[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: [sales@cnbiomedical.com](mailto:sales@cnbiomedical.com)

Mobilephone: 86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.